



1. Pullulan P-800 (MW 853,000)
2. Pullulan P-50 (MW 48,000)
3. Pullulan P-20 (MW 23,700)
4. Pullulan P-5 (MW 5,800)

Column : YMC-Pack Diol-200
 500 × 8.0 mm I.D.
 Eluent : water
 Flow rate : 1.0 mL/min
 Temperature : ambient (24 °C)
 Detection : RI, 16×10^{-6} RIU/FS
 Injection : 180 μ L (0.25 mg/mL)
 Sample : Shodex STANDARD P-82
 manufactured by SHOWA DENKO K.K.