



HPLC - Detector Lamps

Can't find the lamp for your detector or spectrophotometer?

For further details of lamps for UV-Vis Spectrophotometers and other HPLC manufacturers contact us for a complete listing.

PART NUMBER	DESCRIPTION	MODEL	Agilent Technologies (Hewlett-Packard)	MANUFACTURER'S P/N
DHP-901	Deuterium	HP 1040/HP 1050 (G1306A) DAD/HP 1050 DA (1050 MWD)/ HP 1050 MW (79854A) /HP 1090 (75880A) DAD		79883-60002
DHP-902	Deuterium	HP 1080/HP 1081/HP1081B/HP1082B/HP1084/HP1084B		79875-60012
DHP-903	Deuterium	HP 1050 VW (79853C)		79853-60002
DHP-906	Xenon	HP 1046/HP1046A		
DHP-909	Deuterium	HP 8450/8450A		08450-60106
DHP-910	Deuterium	HP 1100 (G1314) VW		G1314-60100
DHP-910LL	LL deuterium	Agilent 1100 VWD long life		
DHP-911	Deuterium	HP 1100 (G1315A) DAD		2140-0590
DHP-911LL	LL deuterium	Agilent 1100 DAD long life		5181-1530
DHP-912	Deuterium	HP 8453		2140-8605
DHP-913	Deuterium	HP 8452 A DAD/HP 8452A Opt 002		08452-60104

PART NUMBER	DESCRIPTION	MODEL	Beckman	MANUFACTURER'S P/N
DAB-905	Deuterium	155		22887153
DAB-906	Deuterium	164/165/167		236920
DAB-909	Deuterium	163		22947029
DAB-913A	Deuterium	Beckman 166 D2 lamp (prealigned) System Gold		538706
DAB-915A	Deuterium	Beckman 168 D2 lamp (prealigned) System Gold		538711
DAB-916	Deuterium	Beckman P/ACE 2000 D2 (prealigned)		

Lamps for Beckman DU Series Spectrophotometers are also available

PART NUMBER	DESCRIPTION	MODEL	Gilson/Rainin	MANUFACTURER'S P/N
DGI-901	Deuterium	Holochrome		HMD(9510)
DGI-902	Deuterium	Spectrochrome		2900-0496
DGI-903	Mercury	111B/112 FW		03-03277-10
DGI-915	Deuterium	115/116/117		100-123
DGI-918	Deuterium	118/119/151/152/155/156		100-326
DGI-918LL	LL deuterium	Gilson 115/116/117/118/119/151 152/155/156 long life		
DGI-919LL	LL deuterium	Gilson 170 DAD D2 long life		
DHP-911	Deuterium	170 Diode Array		2140-60001

PART NUMBER	DESCRIPTION	MODEL	Kontron: Biotek-Kontron	MANUFACTURER'S P/N
DKO-901	Deuterium	LC 430		93-00636
DKO-901H	Deuterium	LC 330		54-02003
DKO-903	Deuterium	Uvikon 720/720LC		035-1001
DKO-904	Xenon	Kontron SFM25 Xe lamp		54-01002
DKO-905	Deuterium	LC 433		93-00636
DKO-906	Deuterium	LC 710/LC 715/Uvikon 820/930		54-02002
DKO-908	Deuterium	Uvikon 725		
DKO-909	Deuterium	Uvikon 740/800/810/860		54-02002
DKO-910	Deuterium	440 DAD		91-910-95
DKO-911	Deuterium	735LC/Uvikon 722/730/922		54-02001
DKO-912	Deuterium	LC 332/335		93-00636
DKO-912H	Deuterium	535 VWD		91-91494

HPLC - Detector Lamps

PART NUMBER	DESCRIPTION	MODEL	Kontron: Biotek-Kontron	MANUFACTURER'S P/N
DKO-912LL	LL deuterium	Kontron 332 335 430 432 433 535 long life		91-91494
DKO-913	Deuterium	Uvikon 941/942/943		54-02002
DKO-913LL	LL deuterium	Kontron 535DAD long life		91-91494
DKO-914LL	LL deuterium	Kontron 540DAD 540+ 545V long life		54-02007
DKO-914	Deuterium	540 DAD/540+ DAD/545DAD		54-02007

PART NUMBER	DESCRIPTION	MODEL	Merck-Hitachi (Hitachi)	MANUFACTURER'S P/N
DHI-901	Deuterium	101/102/111		
DHI-902	Deuterium	100-10/100-40/100-50/100-60		799-9991
DHI-903	Deuterium	150-20/200/220/300/330/340/2000/3000/ 4000/L2500/L3000/L4000/L-4500		885-3570
DHI-908	Deuterium	L4200/L4250/L4500		371140005
DHI-910	Deuterium	LaChrom L4720/L4520/L7400/L450		371140005
DHI-911	Xenon	Hitachi fluorescence detectors F1000/2000/4000 Series		

Lamps for Hitachi Spectrophotometers are also available.

PART NUMBER	DESCRIPTION	MODEL	ABI, Kratos	MANUFACTURER'S P/N
DKA-902	Deuterium	757/769/770/773/775/783/873/FS970/120/130		2900-0496 (-0487) (-0489)
DKA-904	Tungsten	757/769/770/773/775/783/873/FS970/120/130		2450-0213
DKA-910	Deuterium	785A/FS980		2900-0484

PART NUMBER	DESCRIPTION	MODEL	PerkinElmer	MANUFACTURER'S P/N
DPE-903	Deuterium	Lambda 3/7/9		C055-0505
DPE-906	Deuterium	360/460/560		0057-0194
DPE-908	Tungsten	Lambda 2/2S/10/11 and others		
DPE-911	Deuterium	Integral 2000/Integral 4000/LC55/LC65/LC85/LC95		B016-0917/0271-340/0271-1706
DPE-913	Deuterium	LC-90/LC-290		0271-2224
DPE-914	Deuterium	Lambda 2/2S/10/11 and others		B016-0817
DPE-915	Deuterium	Series 200 DAD		N2922010

Lamps for other PerkinElmer Lambda Series Spectrophotometers are also available.

PART NUMBER	DESCRIPTION	MODEL	Shimadzu	MANUFACTURER'S P/N
DSH-901	Deuterium	UV120/UV160/UV160A/UV240/UV260/UV265		200-75503-01
DSH-902	Deuterium	SPD-2A/SPD-3/SPD-4		062-65056
DSH-903	Deuterium	D300L/UV200S		062-70606-00
DSH-912	Xenon	Shimadzu RF530/RF510		
DSH-913	Xenon	Shimadzu RF540/RF535/RF551/RF500		
DSH-914	Xenon	Shimadzu RF1501.5301/5000		
DSH-915	Xenon	RF10A RF10AX		
DSH-916	Deuterium	SPD 6A/SPD-6AV		
DSH-917	Deuterium	SPD 10A/SPD 10AS/SPD-10AV/SPD-10AVP		228-34016-02
DSH-918	Deuterium	SPD-M10AVP PDA		228-34016
DSH-918LL	LL deuterium	Shimadzu SPD-10 Series long life		

HPLC - Detector Lamps


PART NUMBER	DESCRIPTION	MODEL	Thermo Scientific	MANUFACTURER'S P/N
			SpectraPhysics/LDC/Milton Roy/ Linear/Spectronics/TSP	
DLD-901	Deuterium	LDC Spectromonitor III		108035
DLD-907	Deuterium	LDC Spectromonitor I/II		
DLD-909	Deuterium	LDC SM4000 series		900918001
DLD-910	Deuterium	LDC SM3000/3100/3200		108035
DLD-911	Xenon	LDC Fluoromonitor FM4100		
DPY-915	Deuterium	Aquamate/UV1/UV2/UV3/UV4		
DPY-916	Tungsten	Helios Epsilon		
DTL-901	Tungsten	Helios Delta/Gamma		
DSP-901	Deuterium	SP8400/SP8430/SP8440/SP8450/SP8480/SP8490		3302-9540
DSP-903	Deuterium	SP8200/LC871		
DSP-905	Xenon	FL2000/LC304		
DSP-907	Deuterium	SP8480XR/SP8773XR		
DSP-908	Deuterium	Linear UV100/UV200/UV1000/UV2000/ UV3000/Focus/Spectrochrom		9551-0202
DSP-908	Deuterium	Thermo Surveyor		
DSP-912	Deuterium	TSP UV6000LP		108052

PART NUMBER	DESCRIPTION	MODEL	Unicam (Philips/ATI/Pye)	MANUFACTURER'S P/N
DPY-901	Deuterium	4020/4025		4013164-45861
DPY-902	Deuterium	4110		4013166-66428
DPY-903	Deuterium	SP8200/SP8400/SP8600		4013163-75402
DPY-904	Deuterium	600UV		4013161-05993
DPY-905	Deuterium	4021		4013164-39401
DPY-906	Deuterium	8500/8740/8625/8675		9423186-0811
DPY-907	Deuterium	SP8700/SP8750		9423185-03021
DPY-914	Deuterium	4225		9435242-25401
DPY-915	Deuterium	Helios/UV300/UV550		
DTL-901	Tungsten	Helios/UV300/UV550		943-UV9-0001E
DTL-902	Tungsten	SP8200/SP8400/SP8600 (Vis)		9423185

Lamps for Unicam spectrophotometers are also available.

PART NUMBER	DESCRIPTION	MODEL	Varian	MANUFACTURER'S P/N
DVA-901	Deuterium	UV 2050		03-951303-99
DVA-903	Deuterium	UV 50/Varichrom		03-951305
DVA-904	Deuterium	UV100/UV200		03-916077
DVA-905	Deuterium	UV5/2550		00-996800
DVA-906	Deuterium	LC5000/LC5500		03-916077-00
DVA-907	Deuterium	Star 9050		03-916077-00
DVA-909	Deuterium	ProStar 340/345 UV/Vis		03-91615691

Lamps for Varian CARY Spectrophotometers and other models also available.



HPLC - Detector Lamps

PART NUMBER	DESCRIPTION	MODEL	Waters	MANUFACTURER'S P/N
DWA-901	Mercury	440/441/490		WAT 097323
DWA-910	Deuterium	480/481/480LC/481LC/Lambda Max/LC1		WAT 099499
DWA-911	Tungsten	RI/R401/R403/R404		WAT 048419
DWA-912	Cadmium	440/441/490		WAT 097731
DWA-913	Zinc	440/441/490		WAT 097723
DWA-915	Deuterium	484		WAT 080357
DWA-918	Deuterium	486		WAT 080678
DWA-918LC	Deuterium	2486		WAT 080678
DWA-921	Deuterium	996 PDA/2996		WAT 057760
DWA-921LL	LL deuterium	Waters 996		WAT 057760
DWA-923	Xenon	470/475/2475 lamp only		WAT 047018
DWA-926	Deuterium	990/991/994 PDA		WAT 021516
DWA-929	Xenon	474		WAT 047447
DWA-930	Deuterium	2487 Dual Wavelength/2488		WAT 081142
DWA-930LL	LL deuterium	Waters Alliance 2487/2488		WAT 081142

PART NUMBER	DESCRIPTION	MODEL	Cecil	MANUFACTURER'S P/N
DCE-901	Deuterium	Cecil Series 1		2202-0142
DCE-902	Deuterium	Cecil Series 2		2900-0484

Lamps for Cecil UV Spectrophotometers also available

PART NUMBER	DESCRIPTION	MODEL	Dionex (Gynkotek-Softron)	MANUFACTURER'S P/N
DGP-903	Deuterium	Gynkotek UVD-340S		
DDX-901	Deuterium	Dionex Summit PDA-100		

PART NUMBER	DESCRIPTION	MODEL	Knauer	MANUFACTURER'S P/N
DKN-904	Deuterium	8600/8700		
DKN-905	Deuterium	Wellchrom K-2500/2501/2600		

PART NUMBER	DESCRIPTION	MODEL	Severn Analytical	MANUFACTURER'S P/N
DSA-901	Deuterium	SA6500/SA6503/SA6504/SA6508		

ADDITIONAL NOTES

- Deuterium lamps are usually supplied ready for installation with pre-aligned bases. In some cases the holder or bracket from the lamp assembly will need to be reused. Deuterium lamps are usually warranted for 1000 hours of operation. The long-life deuterium lamps are now supplied with an extended operational life of 2000 hours+.
- Xenon lamps are supplied without mounting brackets and require manual alignment.
- Tungsten lamps are pre-aligned.

Agilent 1100
DHP-911LL



Shimadzu
DSH-917LL



Gilson
DGI-918



Varian
DVA-907

